

GOAT ANTI-OAS3 ANTIBODY

SKU: EB12375



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	p100OAS p100 OAS 2-5A synthetase 3 2-5A synthase 3 2'-5'-oligoadenylate synthetase p100 2'-5'-oligoadenylate synthase 3 (2-5')oligo(A) synthetase 3 (2-5')oligo(A) synthase 3 p100 2'-5'-oligoadenylate synthetase 3, 100kDa OAS3
Accession ID	NP_006178.2
Blocking Peptide	EBP12375
Immunogen	Peptide with sequence C-HPIQLDPNQKTPENS, from the internal region of the protein sequence according to NP_006178.2.
Peptide Sequence	C-HPIQLDPNQKTPENS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	4940
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 110kDa band observed in lysates of cell line A431 (calculated MW of 121kDa according to NP_006178.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

