


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-OAZ1 (AA120-132) ANTIBODY





SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias antizyme 1|ODC-Az|OAZ|AZI|ornithine decarboxylase antizyme 1|OAZ1

Names

Immunofluorescence: Strong expression of the protein seen in A431 and U2OS cells.

Usage Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of A431 Summary

cells. Recommended concentration: 10ug/ml.

Accession

NP_004143.1; NP_001287949.1

Blocking

EBP12910 Peptide





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Peptide with sequence C-SRLTDAKRINWRT, from the internal region of the protein sequence according to Immunogen NP_004143.1; NP_001287949.1.

Peptide

C-SRLTDAKRINWRT Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Refrigerated Instructions

Predicted Species

Human

Reactive

Human **Species**

Human

4946 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Western Approx. 25kDa band observed in Human Kidney lysates (calculated MW of 25.4kDa according to NP_004143.1).

Blot Recommended concentration: 0.3-1µg/ml. Primary incubation 1 hour at room temperature.

Application


Pep-ELISA, WB, IF, FC Type


GALLERY IMAGES



Telephone: (650) 697-3600







250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600



