

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-OAZ1 ANTIBODY

SKU: EB12909

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias antizyme 1|ODC-Az|OAZ|AZI|ornithine decarboxylase antizyme 1|OAZ1

Names

Accession NP 004143.1 ID

Blocking

EBP12909 Peptide

Immunogen NP_004143.1. Peptide with sequence C-NSQRDHNLSAN, from the internal region of the protein sequence according to

Peptide C-NSQRDHNLSAN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Dog, Pig, Cow **Species**

Human

4946 **Gene ID**

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Application Pep-ELISA **Type**