

GOAT ANTI-OAZ1 ANTIBODY

SKU: EB12909



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. |
| Unit Size | 100 µg |
| Storage Instructions | Aliquot and store at -20°C. Minimize freezing and thawing. |
| Synonym / Alias Names | antizyme 1 ODC-Az OAZ AZI ornithine decarboxylase antizyme 1 OAZ1 |
| Accession ID | NP_004143.1 |
| Blocking Peptide | EBP12909 |
| Immunogen | Peptide with sequence C-NSQRDHNLSAN, from the internal region of the protein sequence according to NP_004143.1. |
| Peptide Sequence | C-NSQRDHNLSAN |
| Purification Method | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Shipping Instructions | Refrigerated |
| Predicted Species | Human, Dog, Pig, Cow |
| Human Gene ID | 4946 |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png |
| ELISA Detection Limit | Antibody detection limit dilution 1:128000. |
| Application Type | Pep-ELISA |

DOCUMENTS

- [Data Sheet](#)