

GOAT ANTI-OCT6 / SCIP / POU3F1 ANTIBODY

SKU: EB06207



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	octamer-binding transcription factor 6 POU class 3 homeobox 1 SCIP POU domain, class 3, transcription factor 1 OCT-6 OTF6 OCT6 POU3F1
Accession ID	NP_002690.3
Blocking Peptide	EBP06207
Immunogen	Peptide with sequence ATTAQYLPRGPGG-C, from the N Terminus of the protein sequence according to NP_002690.3.
Peptide Sequence	ATTAQYLPRGPGG-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow
Human Gene ID	5453
Mouse Gene ID	18991
Rat Gene ID	192110
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA