



Telephone: (650) 697-3600

## **GOAT ANTI-OCT6 / SCIP / POU3F1 ANTIBODY**

**SKU:** EB06207

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

octamer-binding transcription factor 6|POU class 3 homeobox 1|SCIP|POU domain, class 3, transcription factor **Alias** 

1|OCT-6|OTF6|OCT6|POU3F1 **Names** 

Accession

NP 002690.3 ID

**Blocking** 

EBP06207 **Peptide** 

Peptide with sequence ATTAQYLPRGPGG-C, from the N Terminus of the protein sequence according to Immunogen

NP 002690.3.

**Peptide** 

ATTAQYLPRGPGG-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Cow **Species** 

Human

5453 Gene ID

Mouse **Gene ID** 

18991

**Rat Gene ID** 192110

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:8000.

Limit

**Application** Pep-ELISA

**Type**