


GOAT ANTI-OIP106 / TRAK1 ANTIBODY

SKU: EB06605



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	trafficking protein, kinesin binding 1 OGT(O Glc NAc transferase) interacting protein 106 Kda OGT(O-Glc-NAc transferase)-interacting protein 106 KDa OIP106
Accession ID	NP_001036111.1
Blocking Peptide	EBP06605
Immunogen	Peptide with sequence CGAKLSKQTSRLR, from the C Terminus of the protein sequence according to NP_001036111.1.
Product Comments	This antibody is only expected to recognise isoform 1 (NP_001036111.1) only.
Peptide Sequence	CGAKLSKQTSRLR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	22906
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 90kDa band observed in Human Heart and Human Lung lysates (calculated MW of 77.3kDa according to NP_055780). Recommended concentration: 0.3-1.5µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

