Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-OPRD1 ANTIBODY**

**SKU:** EB09027

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias OPRD|opioid receptor, delta 1|OPRD1

**Names** 

Accession NP 000902.3

**Blocking** 

EBP09027 **Peptide** 

Peptide with sequence TARERVTACTPSDG, from the C Terminus of the protein sequence according to Immunogen

NP 000902.3.

**Peptide** 

TARERVTACTPSDG Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat **Species** 

Human

4985 Gene ID

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:8000.

Limit

Preliminary experiments in Human Brain (Hippocampus, Olfactory bulb) lysates gave no specific signal but low Western background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field -

Blot have any results been reported with other antibodies/lysates?

Application

Pep-ELISA **Type**