

GOAT ANTI-OR7D4 ANTIBODY

SKU: EB08836



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	OR7D4 olfactory receptor, family 7, subfamily D, member 4 OR19-7 OR19-B OR7D4P hg105 odorant receptor family 7 subfamily D member 4 RT olfactory receptor OR19-7 olfactory receptor, family 7, subfamily D, member 4 pseudogene
Accession ID	NP_001005191.1
Blocking Peptide	EBP08836
Immunogen	Peptide with sequence C-SIQARSKDISY , from the internal region of the protein sequence according to NP_001005191.1.
Peptide Sequence	C-SIQARSKDISY
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	125958
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)