



Telephone: (650) 697-3600

GOAT ANTI-ORC4L ANTIBODY

SKU: EB05067





Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

origin recognition complex, subunit 4 homolog|origin recognition complex subunit 4|ORC4P|ORC4|Origin recognition complex, subunit 4, S. cerevisiae, homolog-like|origin recognition complex, subunit 4 (yeast

homolog)-like|HSORC4|Orc4p|origin recognition complex, subunit 4-like (yeast)|ORC4L Names

Accession

NP_002543.2; NP_001177810.1; NP_001177811.1 ID

Blocking

EBP05067 **Peptide**

Immunogen

Peptide with sequence C-DVRQWATSSLSWL, from the C Terminus of the protein sequence according to

NP_002543.2; NP_001177810.1; NP_001177811.1.

Product

This antibody is expected to recognize all reported isoforms (NP_002543.2; NP_001177810.1; NP_001177811.1). Reported variants represent identical protein: NP_002543.2, NP_001177808.1,

Comments

NP_859526.1, NP_859525.1

Peptide

Sequence

C-DVRQWATSSLSWL

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Cow

Species Reactive

Species

Human

Human

5000

Gene ID

Mouse

26428

Gene ID **Rat Gene ID** 295596

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:8000.

Limit

Western

Approx 48kDa band observed in nuclear lysates of cell line Jurkat (calculated MW of 50.4kDa according to Human NP_002543.2 Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour. Preliminary

testing was unsuccessful on cell line NIH-3T3 and Mouse Heart for this particular batch.

Application

Pep-ELISA, WB

Type

Blot





Telephone: (650) 697-3600

GALLERY IMAGES



Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa