



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-OSBP ANTIBODY**

**SKU:** EB06707



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa 50kDa

37kDa

25kDa

20kDa

15kDa

10kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias OSBP1|HGNC:8503|oxysterol binding protein|OSBP

Names

**Accession** 

NP\_002547.1 ID

**Blocking** 

EBP06707 **Peptide** 

Peptide with sequence C-KHQLEETKKEKRTRIP, from the internal region of the protein sequence according to **Immunogen** 

NP 002547.1.

**Peptide** 

C-KHQLEETKKEKRTRIP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Dog, Cow **Species** 

Reactive

Human **Species** 

Human

5007 **Gene ID** 

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:32000.

Limit

Western Approx 85kDa band observed in Human Placenta lysates (calculated MW of 89.4kDa according to NP\_002547).

**Blot** Recommended concentration: 0.1-0.3µg/ml.

**Application** 

Pep-ELISA, WB Type

## **GALLERY IMAGES**



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa