

GOAT ANTI-OSBP ANTIBODY

SKU: EB06707

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	OSBP1 HGNC:8503 oxysterol binding protein OSBP
Accession ID	NP_002547.1
Blocking Peptide	EBP06707
Immunogen	Peptide with sequence C-KHQLEETKKEKRTRIP, from the internal region of the protein sequence according to NP_002547.1.
Peptide Sequence	C-KHQLEETKKEKRTRIP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	5007
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 85kDa band observed in Human Placenta lysates (calculated MW of 89.4kDa according to NP_002547). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

