

GOAT ANTI-OSBP1 (C TERMINUS) ANTIBODY

SKU: EB06149



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	oxysterol binding protein OSBP1 OSBP
Accession ID	NP_002547.1
Blocking Peptide	EBP06149
Immunogen	Peptide with sequence CKEKQDWSSCPDIF, from the C Terminus of the protein sequence according to NP_002547.1.
Peptide Sequence	CKEKQDWSSCPDIF
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	5007
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 100kDa band observed in Human Liver lysates (calculated MW of 89.4kDa according to NP_002547.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

