

GOAT ANTI-P8 / NUPR1 ANTIBODY

SKU: EB06284



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	candidate of metastasis 1 nuclear protein, transcriptional regulator, 1 p8 protein p8 protein (candidate of metastasis 1) COM1 P8 Nuclear Protein 1 NUPR1
Accession ID	NP_001035948.1; NP_036517
Blocking Peptide	EBP06284
Immunogen	Peptide with sequence ATFPPATSAPQQP-C, from the N Terminus of the protein sequence according to NP_001035948.1; NP_036517.
Peptide Sequence	ATFPPATSAPQQP-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	26471
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)