Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-P8 / NUPR1 ANTIBODY**

**SKU:** EB06284

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

candidate of metastasis 1|nuclear protein, transcriptional regulator, 1|p8 protein|p8 protein (candidate of **Alias** 

metastasis 1)|COM1|P8|Nuclear Protein 1|NUPR1 **Names** 

Accession

NP 001035948.1; NP 036517 ID

**Blocking** 

EBP06284 **Peptide** 

Peptide with sequence ATFPPATSAPQQP-C, from the N Terminus of the protein sequence according to Immunogen

NP 001035948.1; NP 036517.

**Peptide** 

ATFPPATSAPQQP-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

**Species** 

Human

Human

26471 Gene ID

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:32000.

Limit

**Application** Pep-ELISA **Type**