

GOAT ANTI-PAIP2 ANTIBODY

SKU: EB05531



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	polyadenylate binding protein-interacting protein PABP-interacting protein 2 PAIP2
Accession ID	NP_057564
Blocking Peptide	EBP05531
Immunogen	Peptide with sequence KDPSRSSTSPSIIN-C, from the N Terminus of the protein sequence according to NP_057564.
Peptide Sequence	KDPSRSSTSPSIIN-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	51247
Mouse Gene ID	67869
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 15kDa band seen in 293 lysate. Recommended for use at 1-3µg/ml.