

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PAIP2 ANTIBODY

SKU: EB05531

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias polyadenylate binding protein-interacting protein|PABP-interacting protein 2|PAIP2

Names

Accession NP_057564 ID

Blocking

EBP05531

Peptide

Peptide with sequence KDPSRSSTSPSIIN-C, from the N Terminus of the protein sequence according to **Immunogen**

NP 057564.

Peptide

KDPSRSSTSPSIIN-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Human, Mouse, Rat **Species**

Reactive

Human **Species**

Human

51247

Gene ID Mouse

67869 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png **Grade**

Western

Western blot: Approx 15kDa band seen in 293 lysate. Recommended for use at 1-3μg/ml. **Blot**