



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PAK1 ANTIBODY

SKU: EB07367








250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

p21/Cdc42/Rac1-activated kinase 1 (yeast Ste20-related)|p21-activated kinase

Alias Names

1|PAKalpha|MGC130001|MGC130000|p21/Cdc42/Rac1-activated kinase 1 (STE20 homolog, yeast)|PAK1

Accession

ID

NP_002567.3

Blocking

Peptide

EBP07367

Peptide with sequence C-NTEKQKKKPKMSDE, from the internal region of the protein sequence according to **Immunogen**

NP 002567.3.

Product Comments

No cross-reactivity exptected with PAK3

Peptide

Sequence

C-NTEKQKKKPKMSDE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

5058

Predicted

Species

Human, Mouse, Rat, Dog

Human, Mouse, Rat

Reactive

Species Human

Gene ID

Mouse 18479 Gene ID

Rat Gene ID 29431

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:1000.

Limit

Western Approx 65kDa band observed in Human Brain (Cerebral Cortex and Hippocampus) lysates (calculated MW of

60.6kDa according to NP_002567.3). Recommended concentration: 0.5-1.5μg/ml. **Blot**

Application

Type

Pep-ELISA, WB



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa