

GOAT ANTI-PALLADIN ANTIBODY

SKU: EB08623



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PALDD palladin, cytoskeletal associated protein CGI-151 FLJ22190 FLJ38193 FLJ39139 KIAA0992 PNCA1 SIH002 palladin sarcoma antigen NY-SAR-77
Accession ID	NP_057165.3
Blocking Peptide	EBP08623
Immunogen	Peptide with sequence SLSDMQEESKNTD-C, from the N Terminus of the protein sequence according to NP_057165.3.
Peptide Sequence	SLSDMQEESKNTD-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	23022
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Preliminary experiments in Human Liver, Heart and Ovary lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)