

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-PAPD5 (AA497-510) ANTIBODY**

**SKU:** EB10938



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

**Alias** 

topoisomerase-related function protein 4-2| terminal uridylyltransferase 3| TUTase 3| PAP-associated domain-

containing protein 5| TRF4-2| FLJ40270| PAP associated domain containing 5|PAPD5

**Accession** 

ID

Names

NP\_001035374.2; NP\_001035375.2

**Blocking** 

**Peptide** 

EBP10938

Peptide with sequence C-TDEVATYRDWISKQ, from the internal region of the protein sequence according to **Immunogen** 

NP\_001035374.2; NP\_001035375.2.

**Product Comments** 

This antibody is expected to recognize both reported isoforms (NP 001035374.2; NP 001035375.2).

**Peptide** 

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

C-TDEVATYRDWISKQ

Shipping

Instructions

Refrigerated

Mouse, Rat

**Predicted** 

**Species** 

Human, Mouse, Rat, Cow

Reactive

**Species** 

Human

64282 Gene ID

Mouse

214627 Gene ID

Rat Gene ID 307745

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

Grade **ELISA** 

**Detection** 

Antibody detection limit dilution 1:32000.

Limit

Western Approx 70kDa band observed in Mouse and Rat Spleen lysates (calculated MW of 69.4kDa according to Mouse

NP\_001157971.1. Recommended concentration: 0.2-0.6µg/ml. **Blot** 

**Application** 

**Type** 

Pep-ELISA, WB





Telephone: (650) 697-3600

## **GALLERY IMAGES**

e erest

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa