

GOAT ANTI-PAPD5 (AA497-510) ANTIBODY

SKU: EB10938



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	topoisomerase-related function protein 4-2 terminal uridylyltransferase 3 TUTase 3 PAP-associated domain-containing protein 5 TRF4-2 FLJ40270 PAP associated domain containing 5 PAPD5
Accession ID	NP_001035374.2; NP_001035375.2
Blocking Peptide	EBP10938
Immunogen	Peptide with sequence C-TDEVATYRDWISKQ, from the internal region of the protein sequence according to NP_001035374.2; NP_001035375.2.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_001035374.2; NP_001035375.2).
Peptide Sequence	C-TDEVATYRDWISKQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow
Reactive Species	Mouse, Rat
Human Gene ID	64282
Mouse Gene ID	214627
Rat Gene ID	307745
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 70kDa band observed in Mouse and Rat Spleen lysates (calculated MW of 69.4kDa according to Mouse NP_001157971.1. Recommended concentration: 0.2-0.6µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa