


GOAT ANTI-PAPIN / PDZK3 ANTIBODY

SKU: EB06236



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names plakophilin-related armadillo repeat protein-interacting PDZ protein-like|activated in prostate cancer protein|PIN1|AIPC|PDZ domain containing 2|likely ortholog of rat plakophilin-related armadillo repeat protein-interacting PDZ protein|PDZ domain containing 2|PDZ domain containing 3|KIAA0300|PAPIN|PDZK3|PDZD2

Usage Summary **Immunofluorescence:** Expression of the protein seen in the plasma membrane and cytoplasm of U2OS cells. Recommended concentration: 10µg/ml.

Accession ID NP_835260.2

Blocking Peptide EBP06236

Immunogen Peptide with sequence PITQDNAVHLHLPLC, from the N Terminus of the protein sequence according to NP_835260.2.

Peptide Sequence PITQDNAVHLHLPLC

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	23037
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Tonsil. Recommended concentration: 6µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IF, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

