


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PAPIN / PDZK3 ANTIBODY

SKU: EB06236



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Alias

Synonym / plakophilin-related armadillo repeat protein-interacting PDZ protein-like|activated in prostate cancer protein|PIN1|AIPC|PDZ domain containing 2|likely ortholog of rat plakophilin-related armadillo repeat proteininteracting PDZ protein|PDZ domain containing 2|PDZ domain containing 3|KIAA0300|PAPIN|PDZK3|PDZD2 Immunofluorescence: Expression of the protein seen in the plasma membrane and

Usage **Summary**

Names

cytoplasm of U2OS cells. Recommended concentration: 10µg/ml.

Accession

ID

NP_835260.2

Blocking Peptide

EBP06236

Immunogen

Peptide with sequence PITQDNAVLHLPLC, from the N Terminus of the protein sequence according to NP 835260.2.

Peptide

Sequence

PITQDNAVLHLPLC

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.





Telephone: (650) 697-3600



Shipping Refrigerated

Instructions

Predicted Species

Human

Reactive

Human **Species**

Human **Gene ID** 23037

Product Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$

IHC Results Paraffin embedded Human Tonsil. Recommended concentration: 6µg/ml.

ELISA

Antibody detection limit dilution 1:32000. Detection

Limit

Type

Application Pep-ELISA, IF, IHC

GALLERY IMAGES

