

GOAT ANTI-PAX6 (C TERMINUS) ANTIBODY

SKU: EB07165



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	paired box homeotic gene-6 paired box homeotic gene 6 (aniridia, keratitis) paired box gene 6 WAGR MGDA MGC17209 D11S812E AN2 AN paired box gene 6 (aniridia, keratitis) PAX6
Accession ID	NP_000271.1; NP_001595.2
Blocking Peptide	EBP07165
Immunogen	Peptide with sequence C-EPDMSQYWRLQ, from the C Terminus of the protein sequence according to NP_000271.1; NP_001595.2.
Product Comments	This antibody is expected to recognise both reported isoforms (NP_000271.1 and NP_001595.2).
Peptide Sequence	C-EPDMSQYWRLQ
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	5080
Mouse Gene ID	18508
Rat Gene ID	25509
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments in human and mouse brain lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA