

# GOAT ANTI-PCBP3 (AA14-28) ANTIBODY

SKU: EB11296



## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias Names</b>	PCBP3  poly(rC)-binding protein 3  poly(rC) binding protein 3  OTTHUMP00000115469  OTTHUMP00000115468  OTTHUMP00000115467  OTTHUMP00000115466  OTTHUMP00000115465  OTTHUMP00000115464  ALPHA-CP3
<b>Accession ID</b>	NP_065389.2; NP_001123613.1
<b>Blocking Peptide</b>	EBP11296
<b>Immunogen</b>	Peptide with sequence PHSTLSTLSHHPQPQ-C, from the N Terminus of the protein sequence according to NP_065389.2; NP_001123613.1.
<b>Product Comments</b>	This antibody is expected to recognize both reported isoforms (NP_065389.2; NP_001123613.1). No cross-reactivity with other PCBP proteins are expected.
<b>Peptide Sequence</b>	PHSTLSTLSHHPQPQ-C
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Human Gene ID</b>	54039
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>ELISA</b>	
<b>Detection Limit</b>	Antibody detection limit dilution 1:32000.
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS

- [Data Sheet](#)