Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-PCBP3 (AA14-28) ANTIBODY**

**SKU:** EB11296

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

**Alias** 

PCBP3| poly(rC)-binding protein 3| poly(rC) binding protein 3| OTTHUMP00000115469| OTTHUMP00000115468| OTTHUMP00000115467| OTTHUMP00000115466| OTTHUMP00000115465| OTTHUMP00000115464| ALPHA-CP3

**Names** 

Accession NP\_065389.2; NP\_001123613.1 ID

**Blocking** 

EBP11296 **Peptide** 

Peptide with sequence PHSTLSTLSHHPQPQ-C, from the N Terminus of the protein sequence according to Immunogen

NP 065389.2; NP 001123613.1.

**Product** This antibody is expected to recognize both reported isoforms (NP\_065389.2; NP\_001123613.1). No cross-

Comments reactivity wth other PCBP proteins are expected.

**Peptide** 

PHSTLSTLSHHPQPQ-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

Method the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

Human

Species Human

54039

**Product** 

**Gene ID** 

Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:32000.

Limit

**Application** Pep-ELISA **Type** 

> Goat Anti-PCBP3 (aa14-28) Antibody https://everestbiotech.com/products/goat-anti-pcbp3-aa14-28-antibody/