



Telephone: (650) 697-3600



GOAT ANTI-PCSK5 ANTIBODY

SKU: EB08888

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

subtilisin/kexin-like protease PC5|protease PC6|proprotein convertase PC5|prohormone convertase **Alias** 5|OTTHUMP00000063532|SPC6|PC6A|PC6|PC5|proprotein convertase subtilisin/kexin type 5|PCSK5 **Names**

Accession

NP 006191.2 ID

Blocking

EBP08888 **Peptide**

Peptide with sequence PDGSYQDTKKNLCRK, from the internal region of the protein sequence according to Immunogen

NP 006191.2.

Peptide

PDGSYQDTKKNLCRK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

Gene ID

5125

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection

Antibody detection limit dilution 1:8000.

Limit

Application Pep-ELISA **Type**