

GOAT ANTI-PCSK5 ANTIBODY

SKU: EB08888



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	subtilisin/kexin-like protease PC5 protease PC6 proprotein convertase PC5 prohormone convertase 5 OTTHUMP0000063532 SPC6 PC6A PC6 PC5 proprotein convertase subtilisin/kexin type 5 PCSK5
Accession ID	NP_006191.2
Blocking Peptide	EBP08888
Immunogen	Peptide with sequence PDGSYQDTKKNLCRK, from the internal region of the protein sequence according to NP_006191.2.
Peptide Sequence	PDGSYQDTKKNLCRK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	5125
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)