


GOAT ANTI-PDE4D3 ANTIBODY

SKU: EB07542



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias phosphodiesterase 4D3|PDE4D3

Names

Usage **Immunofluorescence:** Strong expression of the protein seen in the cytoplasm of A549 cells.

Summary Recommended concentration: 10µg/ml.

Accession ID NP_006194.2

Blocking Peptide EBP07542

Immunogen Peptide with sequence HVNNFPFRRHS-C, from the N Terminus of the protein sequence according to NP_006194.2.

Product Comments	This antibody is expected to recognize isoform 3 (NP_006194.2) only.
Peptide Sequence	HVNNFPFRRHS-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Rat, Pig
Reactive Species	Human
Human Gene ID	5144
Mouse Gene ID	24627
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA, IF

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

