
Telephone: (650) 697-3600








## **GOAT ANTI-PDE4D3 ANTIBODY**

**SKU:** EB07542



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

**Alias** phosphodiesterase 4D3|PDE4D3

**Names** 

Usage <strong>Immunofluorescence:</strong> Strong expression of the protein seen in the cytoplasm of A549 cells. Summary

Recommended concentration: 10µg/ml.

Accession ID

NP\_006194.2

**Blocking Peptide** 

EBP07542

**Immunogen** 

Peptide with sequence HVNNFPFRRHS-C, from the N Terminus of the protein sequence according to

NP\_006194.2.







Telephone: (650) 697-3600

**Product** Comments

This antibody is expected to recognize isoform 3 (NP\_006194.2) only.

**Peptide** 

e erest

**HVNNFPFRRHS-C** 

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.

**Shipping** Instructions

Refrigerated

**Predicted** 

**Species** 

Human, Rat, Pig

Reactive

**Species** 

Human

Human

5144 Gene ID

24627

Mouse

Gene ID

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**Grade ELISA** 

**Detection** 

Antibody detection limit dilution 1:128000.

Limit

**Type** 

**Application** Pep-ELISA, IF

## **GALLERY IMAGES**

