

GOAT ANTI-PDX1 (ZEBRAFISH) ANTIBODY

SKU: EB07974



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	pancreatic and duodenal homeobox gene 1 insulin promoter factor 1, homeodomain transcription factor zgc:100874 pdx-1 ipf1 pancreatic and duodenal homeobox 1 pdx1
Accession ID	NP_571518.2
Blocking Peptide	EBP07974
Immunogen	Peptide with sequence C-QRHPNEDYSQNP, from the internal region of the protein sequence according to NP_571518.2.
Peptide Sequence	C-QRHPNEDYSQNP
Shipping Instructions	Refrigerated
Predicted Species	Zebrafish
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA