
Telephone: (650) 697-3600



e erest

GOAT ANTI-PECAM1 / CD31 (AA218-231) ANTIBODY

SKU: EB12606



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym / platelet endothelial cell adhesion molecule|CD31

antigen|endoCAM|PECAM-1|PECA1|GPIIA'|CD31/EndoCAM|CD31|platelet/endothelial cell adhesion molecule **Alias**

Names 1|PECAM1

Immunofluorescence: Strong expression of the protein seen in the cell junctions/ PM of Usage

Summary Jurkat cells. Recommended concentration: 10µg/ml.

Accession

NP_000433.4 ID

Blocking

EBP12606 **Peptide**

Peptide with sequence C-TSESTKSELVTVTE, from the internal region of the protein sequence according to **Immunogen**

NP 000433.4.

Peptide Sequence

C-TSESTKSELVTVTE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.





Telephone: (650) 697-3600



Shipping Refrigerated Instructions

Predicted

Human

Species Reactive

Human

Species Human

Gene ID

5175

Product Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$

IHC Results Paraffin embedded Human Kidney. Recommended concentration: 4-6µg/ml.

ELISA

Antibody detection limit dilution 1:128000. Detection

Limit

Type

Application Pep-ELISA, IF, IHC

GALLERY IMAGES

