


GOAT ANTI-PECAM1 / CD31 (AA218-231) ANTIBODY

SKU: EB12606



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names platelet endothelial cell adhesion molecule|CD31 antigen|endoCAM|PECAM-1|PECA1|GPIIA'|CD31/EndoCAM|CD31|platelet/endothelial cell adhesion molecule 1|PECAM1

Usage Summary Immunofluorescence: Strong expression of the protein seen in the cell junctions/ PM of Jurkat cells. Recommended concentration: 10µg/ml.

Accession ID NP_000433.4

Blocking Peptide EBP12606

Immunogen Peptide with sequence C-TSESTKSELVTVE, from the internal region of the protein sequence according to NP_000433.4.

Peptide Sequence C-TSESTKSELVTVE

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	5175
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Kidney. Recommended concentration: 4-6µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA, IF, IHC

GALLERY IMAGES

