

GOAT ANTI-PECANEX HOMOLOG ANTIBODY

SKU: EB05009



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FLJ45663 FLJ23409 PCNXL1 pecanex homolog (Drosophila) KIAA0805 pecanex-like 1 KIAA0995 PCNX
Accession ID	NP_055797.2
Blocking Peptide	EBP05009
Immunogen	Peptide with sequence C-TVIETGVLELGAEV, from the C Terminus of the protein sequence according to NP_055797.2.
Peptide Sequence	C-TVIETGVLELGAEV
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow, Pig
Human Gene ID	22990
Mouse Gene ID	54604
Rat Gene ID	314288
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments gave no signal but low background in Human Skeletal Muscle and ovary lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Application
Type Pep-ELISA

DOCUMENTS

- [Data Sheet](#)