Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-PECANEX HOMOLOG ANTIBODY**

**SKU:** EB05009

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias FLJ45663|FLJ23409|PCNXL1|pecanex homolog (Drosophila)|KIAA0805|pecanex-like 1|KIAA0995|PCNX

**Names** 

Accession NP 055797.2 ID

**Blocking** 

EBP05009 **Peptide** 

Peptide with sequence C-TVIETGVLELGAEV, from the C Terminus of the protein sequence according to Immunogen

NP 055797.2.

**Peptide** 

C-TVIETGVLELGAEV Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Dog, Cow, Pig **Species** 

Human

22990 Gene ID

Mouse

54604 **Gene ID** 

Rat Gene ID 314288

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:64000.

Limit

Blot

Preliminary experiments gave no signal but low background in Human Skeletal Muscle and ovary lysates at up to Western

1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

**Application Type**Pep-ELISA