

GOAT ANTI-PEPP2 / PLEKHA5 ANTIBODY

SKU: EB06165



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	phosphoinositol 3-phosphate-binding protein-2 FLJ31492 KIAA1686 FLJ10667 PEPP2 pleckstrin homology domain containing, family A member 5 PLEKHA5
Accession ID	NP_061885.2; NP_001137293.1
Blocking Peptide	EBP06165
Immunogen	Peptide with sequence QPQLTEGSHFMCV, from the C Terminus of the protein sequence according to NP_061885.2; NP_001137293.1.
Peptide Sequence	QPQLTEGSHFMCV
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	54477
Mouse Gene ID	109135
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Preliminary experiments in Human Heart and Kidney lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA