

GOAT ANTI-PEPP3 / PLEKHA6 ANTIBODY

SKU: EB06164



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	phosphoinositol 3-phosphate-binding protein 3 MGC176733 pleckstrin homology domain containing, family A member 6 phosphoinositol 3-phosphate-binding protein-3 KIAA0969 PEPP3 PLEKHA6
Accession ID	NP_055750.2
Blocking Peptide	EBP06164
Immunogen	Peptide with sequence C-SPRGADSSYTMRV, from the C Terminus of the protein sequence according to NP_055750.2.
Peptide Sequence	C-SPRGADSSYTMRV
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	22874
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	No signal obtained yet but a preliminary experiment showed low background in Human Heart extract at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)