Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PEPP3 / PLEKHA6 ANTIBODY

SKU: EB06164

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

phosphoinositol 3-phosphate-binding protein 3|MGC176733|pleckstrin homology domain containing, family A **Alias** member 6|phosphoinositol 3-phosphate-binding protein-3|KIAA0969|PEPP3|PLEKHA6

Names

Accession NP 055750.2

Blocking

EBP06164 **Peptide**

Peptide with sequence C-SPRGADSSYTMRV, from the C Terminus of the protein sequence according to Immunogen

NP 055750.2.

Peptide

C-SPRGADSSYTMRV Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

22874 Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western

No signal obtained yet but a preliminary experiment showed low background in Human Heart extract at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other

Blot antibodies/lysates?

Application

Pep-ELISA **Type**