



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PEPT2 ANTIBODY

SKU: EB09197

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias SLC15A2|solute carrier family 15 (H+/peptide transporter), member 2|PEPT2

Names

Accession NP 066568.2 ID

Blocking

EBP09197 **Peptide**

Peptide with sequence C-QGNMIKLETKKTK, from the C Terminus of the protein sequence according to Immunogen

NP 066568.2.

Peptide

C-QGNMIKLETKKTK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Pig **Species**

Human

6565 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:8000.

Limit

Preliminary experiments in Human Brain (Frontal and Cerebral Cortex) and Kidney lysates gave no specific signal Western but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in

Blot the field - have any results been reported with other antibodies/lysates?

Application

Pep-ELISA **Type**