



GOAT ANTI-PEPT2 ANTIBODY

SKU: EB09197



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SLC15A2 solute carrier family 15 (H+/peptide transporter), member 2 PEPT2
Accession ID	NP_066568.2
Blocking Peptide	EBP09197
Immunogen	Peptide with sequence C-QGNMIKLETKTK, from the C Terminus of the protein sequence according to NP_066568.2.
Peptide Sequence	C-QGNMIKLETKTK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Pig
Human Gene ID	6565
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Preliminary experiments in Human Brain (Frontal and Cerebral Cortex) and Kidney lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

DOCUMENTS

- [Data Sheet](#)