



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PER1 ANTIBODY

SKU: EB06842

SPECIFICATIONS

Unit Size

100 μg

Synonym / Alias

period 1|period (Drosophila) homolog 1|circadian pacemaker protein RIGUI|Period, drosophila, homolog

of|hPER|RIGUI|PER|MGC88021|HGNC:8845|period homolog 1 (Drosophila)|PER1 Names

Accession

ID

NP_002607.2

Blocking Peptide

EBP06842

Peptide with sequence CRSKAKRSRHHQN, from the internal region of the protein sequence according to

Peptide

Sequence

CRSKAKRSRHHQN

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Pig, Cow

Human

Gene ID

5187

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Antibody detection limit dilution 1:16000.

Limit

Western Blot

Preliminary experiments in Human Brain lysates gave no specific signal but low background (at antibody

concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been

reported with other antibodies/lysates?

Application

Type

Pep-ELISA