

GOAT ANTI-PERP ANTIBODY

SKU: EB07234

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	p53-induced protein PIGPC1 keratinocytes associated protein 1 RP3-496H19.1 1110017A08Rik dj496H19.1 THW PIGPC1 KRTCAP1 KCP1 HGNC:17637 PERP, TP53 apoptosis effector PERP
Accession ID	NP_071404.2
Blocking Peptide	EBP07234
Immunogen	Peptide with sequence CLPNYEDDLLGNAK, from the C Terminus of the protein sequence according to NP_071404.2.
Peptide Sequence	CLPNYEDDLLGNAK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	64065
Mouse Gene ID	64058
Rat Gene ID	292949
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Skin. Recommended concentration: 5µg/ml.
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 20kDa band observed in Human Liver lysates (calculated MW of 21.4kDa according to NP_071404.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

