



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-PERP ANTIBODY**

**SKU:** EB07234



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa 20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

p53-induced protein PIGPC1|keratinocytes associated protein

Synonym /

**Alias** 1|RP3-496H19.1|1110017A08Rik|dJ496H19.1|THW|PIGPC1|KRTCAP1|KCP1|HGNC:17637|PERP, TP53 apoptosis

effector|PERP Names

Accession

NP\_071404.2

**Blocking** 

ID

EBP07234 **Peptide** 

Peptide with sequence CLPNYEDDLLGNAK, from the C Terminus of the protein sequence according to **Immunogen** 

NP 071404.2.

**Peptide** 

**CLPNYEDDLLGNAK** Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Dog **Species** 

Reactive

Human **Species** 

Human

64065 Gene ID

Mouse

64058 Gene ID

**Rat Gene ID** 292949

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png **Grade** 

IHC Results Paraffin embedded Human Skin. Recommended concentration: 5µg/ml.

**ELISA** 

**Detection** Antibody detection limit dilution 1:64000.

Limit

Western Approx 20kDa band observed in Human Liver lysates (calculated MW of 21.4kDa according to NP\_071404.2).

**Blot** Recommended concentration: 1-3µg/ml.

**Application** 

Pep-ELISA, WB, IHC **Type** 



Telephone: (650) 697-3600



## **GALLERY IMAGES**

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa 20kDa

15kDa

