


GOAT ANTI-PEX13 ANTIBODY

SKU: EB05899



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias peroxisomal biogenesis factor 13|NALD|ZWS|peroxisome biogenesis factor 13|PEX13

Names

Accession ID NP_002609.1

Blocking Peptide EBP05899

Immunogen Peptide with sequence C-DSIGKDGEKQDL, from the C Terminus of the protein sequence according to NP_002609.1.

Peptide Sequence C-DSIGKDGEKQDL

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human
Reactive Species	Human
Human Gene ID	5194
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Liver shows vesicular staining in the cytoplasm of the hepatocytes. Recommended concentration, 2-4µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

