




Telephone: (650) 697-3600



## **GOAT ANTI-PEX13 ANTIBODY**

**SKU:** EB05899



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

peroxisomal biogenesis factor 13|NALD|ZWS|peroxisome biogenesis factor 13|PEX13 **Alias** 

**Names** 

Accession NP\_002609.1 ID

**Blocking** 

EBP05899 **Peptide** 

**Immunogen** 

Peptide with sequence C-DSIGKDGEKQDL, from the C Terminus of the protein sequence according to

NP\_002609.1.

**Peptide** 

C-DSIGKDGEKQDL

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

**Predicted** 

Human

**Species** 

Reactive **Species** 

Human

Human

Gene ID

5194

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

In paraffin embedded Human Liver shows vesicular staining in the cytoplasm of the hepatocytes. **IHC Results** Recommended concentration, 2-4µg/ml.

**ELISA** 

Detection

Antibody detection limit dilution 1:8000.

Limit

Type

**Application** Pep-ELISA, IHC

## **GALLERY IMAGES**

