



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-PEX26 ANTIBODY**

**SKU:** EB10092







250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

PEX26|Pex26pM1T|PEX26M1T|peroxisome biogenesis fac|peroxisome biogenesis disorder, complementation

**Alias** group A|peroxisome biogenesis disorder, complementation group 8|peroxisome assembly protein

26|peroxisomal biogenesis factor 26|peroxin-26|FLJ20695 Names

**Accession** 

NP\_060399.1

**Blocking** 

EBP10092 **Peptide** 

Peptide with sequence C-QKPNLEGSVSHK, from the internal region of the protein sequence according to **Immunogen** 

NP 060399.1.

**Product** 

Reported variants represent identical protein: NP 060399.1, NP 001121121.1

**Comments Peptide** 

Sequence

C-QKPNLEGSVSHK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive

Human **Species** 

Human

55670

Gene ID

**Product** 

Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:64000.

Limit Western

Approx 33kDa band observed in Human Kidney lysates (calculated MW of 33.9kDa according to NP\_060399.1).

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

Recommended concentration: 0.5-1.5µg/ml. Blot

**Application** 

Type

Pep-ELISA, WB

## **GALLERY IMAGES**

Telephone: (650) 697-3600



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa