



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-PEX6 ANTIBODY**

**SKU:** EB10087



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

PEX6|PXAAA1|peroxisome biogenesis factor 6|peroxisomal biogenesis factor 6|peroxisomal AAA-type ATPase

**Alias** 1|PAF-2|PAF2 Names

**Accession** 

NP\_000278.3; NP\_001303242.1 ID

**Blocking** 

EBP10087 **Peptide** 

Peptide with sequence C-EQELLRYKRIQRK, from the C Terminus of the protein sequence according to **Immunogen** 

NP\_000278.3; NP\_001303242.1.

**Peptide** C-EQELLRYKRIQRK

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Cow **Species** 

Reactive

Rat **Species** 

Human

5190 Gene ID

Mouse

224824

**Gene ID** 

**Rat Gene ID** 117265

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**Grade ELISA** 

**Detection** Antibody detection limit dilution 1:128000.

Limit

Western Approx 100kDa band observed in Rat Liver lysates (calculated MW of 104kDa according to Rat NP 476466.1).

Blot Recommended concentration: 0.3-1µg/ml.

**Application** 

Pep-ELISA, WB

**Type** 



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GALLERY IMAGES**

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa