

GOAT ANTI-PEX6 ANTIBODY

SKU: EB10087



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PEX6 PXAAA1 peroxisome biogenesis factor 6 peroxisomal biogenesis factor 6 peroxisomal AAA-type ATPase 1 PAF-2 PAF2
Accession ID	NP_000278.3; NP_001303242.1
Blocking Peptide	EBP10087
Immunogen	Peptide with sequence C-EQELLRYKRIQRK, from the C Terminus of the protein sequence according to NP_000278.3; NP_001303242.1.
Peptide Sequence	C-EQELLRYKRIQRK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow
Reactive Species	Rat
Human Gene ID	5190
Mouse Gene ID	224824
Rat Gene ID	117265
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 100kDa band observed in Rat Liver lysates (calculated MW of 104kDa according to Rat NP_476466.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

