

GOAT ANTI-PGBD1 ANTIBODY

SKU: EB08242



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	cerebral protein-4 OTTHUMP00000039508 dj874C20.4 SCAND4 HUCEP-4 piggyBac transposable element derived 1 PGBD1
Accession ID	NP_115896.1
Blocking Peptide	EBP08242
Immunogen	Peptide with sequence C-EKWSHLSLTRRN, from the internal region of the protein sequence according to NP_115896.1.
Peptide Sequence	C-EKWSHLSLTRRN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	84547
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)