

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PGBD2 ANTIBODY

SKU: EB08243

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias hypothetical protein LOC267002|piggyBac transposable element derived 2|PGBD2

Names

Accession NP 001017434.1; NP 733843.1 ID

Blocking

EBP08243 **Peptide**

Peptide with sequence C-ESNADTTSQGRRSRR, from the internal region of the protein sequence according to Immunogen

NP_001017434.1; NP_733843.1 .

Product This antibody is expected to recognise both reported isoforms (NP_001017434.1; NP_733843.1) Comments

Peptide

C-ESNADTTSQGRRSRR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human Species

Human

267002

Gene ID Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Application Pep-ELISA **Type**

> Goat Anti-PGBD2 Antibody https://everestbiotech.com/products/goat-anti-pgbd2-antibody/